

ALPHA PLAN



Product brochure
alphaFLOW[®] BASIC

alphaFLOW[®] BASIC

Wet density test

The heart of every filtration product is the filter membrane. Before it can be released for use, it must successfully pass a strictly monitored integrity test – especially in the case of medical devices.

The alphaFLOW[®] Basic guarantees 100 % reliable production testing. In the process, the filter membrane is tested for leakage with high-purity water. The tests are carried out in a membrane-friendly manner and defective products are reliably separated.

Note: The alphaFlow[®] Pro variant combines the processes of the alphaFLOW[®] Basic and the alphaMICRO[®]. In this way, you cover all processes with one machine and do not have to install the alphaMICRO[®] in addition.

	BASIC	PRO
100 % reliable testing for membrane integrity	✓	✓
Integrated automatic plant sanitisation (superheated steam)	✓	✓
BDE system	✓	✓
Fully automatic plant sanitisation	✗	✓
Drying of the products directly after wet testing	✗	✓
Drying of the products directly after wet testing	✗	✓

Technology processes

